
Table of Contents

.....	1
LOAD IN THE *.bdf FILE OF THE CURRENT SUBJECT	1
*****EXTRACT PHOTODIODE	
ONSETS*****	6
FILTER PHOTODIODE SIGNALS TO REMOVE HIGHER FREQUENCY OSCILLATORY	
ACTIVITY.	6
DETREND BY FITTING POLYNOMIAL CURVE.	9
HALF WAVE RECTIFY (HWR) THE FILTERED SIGNAL AND INVERT.	12
FIND ONSETS (DIFF(d_hwr==1) AND OFFSETS (diff(D_hwr==-1).	15
PLOT THE DETRENDED AND FILTERED PHOTODIODE SIGNAL (RIGHT) AND SIGNAL	
CONVERTED TO STEP FUNCTION (LEFT).	16
PLOT THE DETAIL OF THE PHOTODIODE SIGNALS TO SHOW DETECTED ONSETS	
AND OFFSETS.	19
CATEGORIZE THE FEEDBACK ONSET TIMES BASED ON DURATIONS	22
DETECT THE TRIGGERS BASED ON VIDEO-TYPE CORRESPONDING TO THE CUR-	
RENT BLOCK.	25
CALCULATE ITI AND DETERMINE FEEDBACKS FALLING INTO THE SAME TRIAL.	29
ADD TRIAL NUMBER INFORMATION TO THE EVENT FIELD OF THE EEG STRUCTURE	
(EEGLAB based).	30
SAVE ALL THE INFORMATION INTO A MAT FILE AS A TABLE FOR EACH FILM	33

```
% Date: 29-11-2019                Programmed by: D. Bolger
% Script to detect the onset times of the feedbacks of interest and to
% determine the modality and type of each feedback.
% Script applied in the Brain-IHM project.
% *****
```

LOAD IN THE *.bdf FILE OF THE CURRENT SUBJECT

Loads in the *.bdf file, ensuring that 74 channels are included so that the ERGO1 and ERGO2 data are loaded.

```
curr_suj = 's01';
Dirbase =
    fullfile(filesep,'Users','bolger','Documents','work','Projects','Project-
PEPs','PEPs_raw_Subjects',curr_suj,filesep);
Savebase =
    fullfile(filesep,'Users','bolger','Documents','work','Projects','Project-
PEPs','PEPs_Preprocess_Subjects',curr_suj,filesep);
case_spec = 'S01';    % Subject 1 videos 1 and 2 are recorded on same
    *.bdf file.

allfiles= dir(Dirbase);
fileIndex = find(~[allfiles.isdir]);
filenum = dir(strcat(Dirbase,'*.bdf'));                                %find all
    the *.bdf files in the current folder
```

```

filenom = {filenum.name};
ergsig1 = cell(1,length(filenom));
ergsig2 = cell(1,length(filenom));

% OPEN EEGLAB SESSION
[ALLEEG, EEG, CURRENTSET, ALLCOM] = eeglab;
[ALLEEG, EEG] = eeg_store(ALLEEG, EEG, CURRENTSET);
filenoms01 = [1 1 1 1];

for scnt = 1:length(filenoms01)

    % The following three lines is added to resolve a bug occurring
    when
    % opening the *.bdf file.
    x = fileparts( which('sopen') );
    rmpath(x);
    addpath(x, '-begin');

    if strcmp(filenom{1,1}(1:3),case_spec) % Needed to add this to
deal with *.bdf file of two videos.

        if scnt ==1
            fullDir = strcat(Dirbase,filenom{1,scnt});
            fnom = filenom{1,scnt}(1:end-4);
            EEG = pop_biosig(fullDir, 'channels',[1:76],'blockrange',
[0 275], 'ref', [], 'refoptions',{'keepref' 'off'} );
            ergsig1{1,scnt} = EEG.data(75,:);
            ergsig2{1,scnt} = EEG.data(76,:);

            EEG = pop_select( EEG,'nochannel',{'Erg1' 'Erg2'});
            EEG = eeg_checkset( EEG );

            [ALLEEG, EEG, CURRENTSET] = pop_newset(ALLEEG, EEG,
CURRENTSET,'setname',char(fnom),'gui','off'); % Create a new dataset
for the current raw datafile
            [ALLEEG, EEG, CURRENTSET] = eeg_store(ALLEEG, EEG,
CURRENTSET);
            EEG = eeg_checkset( EEG );
            EEG =
pop_saveset( EEG, 'filename',char(fnom),'filepath',Savebase);
            % Saves a copy of the current resampled dataset to the current
directory
            eeglab redraw

        elseif scnt ==2

            counttemp = 1;
            fullDir = strcat(Dirbase,filenom{1,1});
            fnom = 'S01_List1a_Film2';
            EEG = pop_biosig(fullDir, 'channels',[1:76],'blockrange',
[280 477], 'ref', [], 'refoptions',{'keepref' 'off'} );
            ergsig1{1,scnt} = EEG.data(75,:);
            ergsig2{1,scnt} = EEG.data(76,:);

```

```

        EEG = pop_select( EEG, 'nochannel', {'Erg1' 'Erg2'});
        EEG = eeg_checkset( EEG );

        [ALLEEG, EEG, CURRENTSET] = pop_newset(ALLEEG, EEG,
        CURRENTSET, 'setname', char(fnom), 'gui', 'off'); % Create a new dataset
        for the current raw datafile
        [ALLEEG, EEG, CURRENTSET] = eeg_store(ALLEEG, EEG,
        CURRENTSET);
        EEG = eeg_checkset( EEG );
        EEG =
        pop_saveset( EEG, 'filename', char(fnom), 'filepath', Savebase);
        % Saves a copy of the current resampled dataset to the current
        directory
        eeglab redraw

    elseif scnt>2
        counttemp = counttemp+1;
        fullDir = strcat(Dirbase, filenom{1, counttemp});
        fnom = filenom{1, counttemp}(1:end-4);
        EEG = pop_biosig(fullDir, 'channels', [1:76], 'ref',
        [], 'refoptions', {'keepref' 'off'} );

        ergsig1{1, scnt} = EEG.data(75, :);
        ergsig2{1, scnt} = EEG.data(76, :);

        EEG = pop_select( EEG, 'nochannel', {'Erg1' 'Erg2'});
        EEG = eeg_checkset( EEG );

        [ALLEEG, EEG, CURRENTSET] = pop_newset(ALLEEG, EEG,
        CURRENTSET, 'setname', char(fnom), 'gui', 'off'); % Create a new dataset
        for the current raw datafile
        [ALLEEG, EEG, CURRENTSET] = eeg_store(ALLEEG, EEG,
        CURRENTSET);
        EEG = eeg_checkset( EEG );
        EEG =
        pop_saveset( EEG, 'filename', char(fnom), 'filepath', Savebase);
        % Saves a copy of the current resampled dataset to the current
        directory
        eeglab redraw
    end

    else % For all other subjects who have 4 *.bdf files each.

        fullDir = strcat(Dirbase, filenom{1, scnt});
        fnom = filenom{1, scnt}(1:end-4);

        % Opening up *.bdf file and saving as a *.set file.
        EEG = pop_biosig(fullDir, 'channels', [1:76], 'ref',
        [], 'refoptions', {'keepref' 'off'} );

        ergsig1{1, scnt} = EEG.data(75, :);

```

```

    ergsig2{1,scnt} = EEG.data(76,:);

    EEG = pop_select( EEG,'nochannel',{ 'Erg1' 'Erg2' });    %
    The photodiode signals are on ERG1 and ERG2- need to extract these
    channels.
    EEG = eeg_checkset( EEG );

    [ALLEEG, EEG, CURRENTSET] = pop_newset(ALLEEG, EEG,
    CURRENTSET,'setname',char(fnom),'gui','off'); % Create a new dataset
    for the current raw datafile
    [ALLEEG, EEG, CURRENTSET] = eeg_store(ALLEEG, EEG,
    CURRENTSET);
    EEG = eeg_checkset( EEG );
    EEG =
    pop_saveset( EEG, 'filename',char(fnom),'filepath',Savebase);
    % Saves a copy of the current resampled dataset to the current
    directory
    eeglab redraw
    end % end of if

end

eeglab: options file is ~/eeg_options.m
EEGLAB: adding "Biosig" to the path; subfolders (if any) might be
missing from the path
EEGLAB: adding "ERPWAVELABv" to the path; subfolders (if any) might be
missing from the path
EEGLAB: adding "Fileio" to the path; subfolders (if any) might be
missing from the path
EEGLAB: adding "dipfit" v2.3 (see >> help eegplugin_dipfit)
EEGLAB: adding "erplab" v7.0.0 (see >> help eegplugin_erplab)
EEGLAB: adding "firfilt" v1.6.1 (see >> help eegplugin_firfilt)
EEGLAB: adding "fullRankAveRef" v0.10 (see >> help
eegplugin_fullRankAveRef)
EEGLAB: adding "limo_eeg-master" v2.0 (see >> help eegplugin_limo)
EEGLAB: adding "loreta" v1.1 (see >> help eegplugin_loreta)
EEGLAB: adding "rERP" v0.4 (see >> help eegplugin_rerp)
Creating a new ALLEEG dataset 1
Warning: line (423: 65152,"V/s",,"",,#obsolete#) not valid
Warning: line (428: 65312,"mHg s-1",,"",,#obsolete#) not valid
Warning: line (430: 65376,"r.p.m",,"rotations per minute",#obsolete#)
not valid
Warning: line (431: 65408,"B", "Bel", "relative power
decibel",#obsolete#) not valid
Warning: line (432: 65440,"dyne s m2 cm-5","Vascular Resistance
Index","dyne seconds square meter per centimetre to the power of
5",#obsolete#) not valid
Warning: line (433: 65440,"dyne*s*mÂ²/cm^5","Vascular Resistance
Index","dyne seconds square meter per centimetre to the power of
5",#obsolete#) not valid
Warning: line (436: 65504,"T",,"<magnitude> Tesla",#obsolete#) not
valid
Reading data in BDF format...

```

```

eeg_checkset note: upper time limit (xmax) adjusted so (xmax-
xmin)*srate+1 = number of frames
Importing data events...
eeg_checkset note: creating the original event table (EEG.urevent)
eeg_checkset note: re-creating the original event table (EEG.urevent)
Removing 2 channel(s)...
Creating a new ALLEEG dataset 1
Saving dataset...
Reading data in BDF format...
eeg_checkset note: upper time limit (xmax) adjusted so (xmax-
xmin)*srate+1 = number of frames
Importing data events...
Removing 2 channel(s)...
Creating a new ALLEEG dataset 2
Saving dataset...
Reading data in BDF format...
eeg_checkset note: upper time limit (xmax) adjusted so (xmax-
xmin)*srate+1 = number of frames
Importing data events...
eeg_checkset note: creating the original event table (EEG.urevent)
eeg_checkset note: re-creating the original event table (EEG.urevent)
Removing 2 channel(s)...
Creating a new ALLEEG dataset 3
Saving dataset...
Reading data in BDF format...
eeg_checkset note: upper time limit (xmax) adjusted so (xmax-
xmin)*srate+1 = number of frames
Importing data events...
eeg_checkset note: creating the original event table (EEG.urevent)
eeg_checkset note: re-creating the original event table (EEG.urevent)
Removing 2 channel(s)...
Creating a new ALLEEG dataset 4
Saving dataset...

```

#4: S01_List1a_Film4

```

Filename: ...ts/s01/S01 List1a Film4.set
Channels per frame          74
Frames per epoch            477184
Epochs                      1
Events                      1
Sampling rate (Hz)          2048
Epoch start (sec)           0.000
Epoch end (sec)             233.000
Reference                    unknown
Channel locations            No (labels only)
ICA weights                  No
Dataset size (Mb)            145.1

```

*****EXTRACT PHOTODIODE ONSETS*****

Need to carry out some pre-processing on the photodiode signals to facilitate the detection of their onsets. The idea is to convert the photodiode signals to step functions.

```
allfiles2 = dir(Savebase); % Finding the number of
*.set files in the current pre-processed data folder.
filenum2 = dir(strcat(Savebase, '*.set')); % Find all the *.set files
in the current folder
filenom2 = {filenum2.name};
fb_sessions = filenom2; % Contains the name of each
video type presented.

% Set parameters for photodiode signal pre-processing.
thresh_val = 0;
fs = EEG.srate; % Sampling frequency.
[b,a] = butter(2,6./(fs/2)); % Low pass filter with 6Hz cutoff to
deal with photodiode flickering- get coefficients.

% Initialise variables.
onoffsets_ergo1 = cell(1,size(ergsig1,2));
onoffsets_ergo2 = cell(1,size(ergsig1,2));
onsets_ergo1 = cell(1,size(ergsig1,2));
onsets_ergo2 = cell(1,size(ergsig1,2));
offsets_ergo1 = cell(1,size(ergsig1,2));
offsets_ergo2 = cell(1,size(ergsig1,2));
Times_bloc = cell(length(fb_sessions),1); % Time vector for each
video.

for ergcnt = 1:size(ergsig1,2)

    Times_bloc{ergcnt,1} = ALLEEG(ergcnt).times;
    time = Times_bloc{ergcnt,1};
```

FILTER PHOTODIODE SIGNALS TO REMOVE HIGHER FREQUENCY OSCILLATORY ACTIVI- TY.

```
ergsig1_curr = ergsig1{1,ergcnt};
ergsig2_curr = ergsig2{1,ergcnt};
ergsig1D = detrend(ergsig1_curr,0); % Detrend the photodiode1
ergsig2D = detrend(ergsig2_curr,0); % Detrend the photodiode2


ergsig1Dfilt = filtfilt(b,a,double(ergsig1D));
ergsig2Dfilt = filtfilt(b,a,double(ergsig2D));



% Plot the detrended and filtered signals
figure('Name',[fb_sessions{1,ergcnt},': detrended +
filtered'],'NumberTitle','off');
```


```

subplot(1,2,1)
plot(time, ergsig1Dfilt); title('Photodiode1 Signal (ERG1):
Detrend+Filtered');
subplot(1,2,2)
plot(time, ergsig2Dfilt); title('Photodiode2 Signal (ERG2):
Detrend+Filtered');

```







DETREND BY FITTING POLYNOMIAL CURVE.

```


p1 = polyfit(1:length(ergsig1Dfilt),ergsig1Dfilt,1); %
Photodiode1
f1 = polyval(p1,1:length(ergsig1Dfilt));
ergsig1Dfilt_corr = ergsig1Dfilt - f1;

p2 = polyfit(1:length(ergsig2Dfilt),ergsig2Dfilt,1); %
Photodiode2
f2 = polyval(p2,1:length(ergsig2Dfilt));
ergsig2Dfilt_corr = ergsig2Dfilt - f2;


% Plot the detrended using polynomial fitting and filtered signals
figure('Name',[fb_sessions{1,ergcnt},': detrended poly
+filtered'],'NumberTitle','off');
subplot(1,2,1)
plot(time, ergsig1Dfilt_corr); title('Photodiode1 Signal (ERG1):
Detrend poly+Filter');
subplot(1,2,2)
plot(time, ergsig2Dfilt_corr); title('Photodiode2 Signal (ERG2):
Detrend poly+Filter');



```

Photodiode1 Signal (ERG1): Detrend poly+Filter **Photodiode2 Signal (ERG2): Detrend poly+Filter**



Photodiode1 Signal (ERG1): Detrend poly+Filter **Photodiode2 Signal (ERG2): Detrend poly+Filter**





HALF WAVE RECTIFY (HWR) THE FILTERED SIGNAL AND INVERT.





```
ergsig1_hwr = zeros(size(ergsig1Dfilt_corr));
ipos = find(ergsig1Dfilt_corr>0);
ergsig1_hwr(ipos)= ergsig1Dfilt_corr(ipos).*1;





ergsig2_hwr = zeros(size(ergsig2Dfilt_corr));
ipos1 = find(ergsig2Dfilt_corr<0);
ergsig2_hwr(ipos1) = ergsig2Dfilt_corr(ipos1).*-1;

% Set all activity <= 1/3 of max to 0.
ergsig1_hwr(ergsig1_hwr<=max(ergsig1_hwr)/3) = 0; ergsig1_hwrz =
ergsig1_hwr;
ergsig2_hwr(ergsig2_hwr<=max(ergsig2_hwr)/3) = 0; ergsig2_hwrz =
ergsig2_hwr;

% Turn all trigger signals into step functions.
ergsig1_hwrst = ergsig1_hwrz;
ergsig2_hwrst = ergsig2_hwrz;
ergsig1_hwrst( ergsig1_hwrz>0) = 1;
ergsig2_hwrst( ergsig2_hwrz>0) = 1;

% Plot the detrended using polynomial fitting, filtered, hwr and
step-function converted signals
figure('Name',[fb_sessions{1,ergcnt},': detrended, filtered, hwr,
step'],'NumberTitle','off');
subplot(1,2,1)
plot(time, ergsig1_hwrst); title('Photodiode1 Signal (ERG1):
HWR'); set(gca,'YLim',[0 1.5]);
subplot(1,2,2)
plot(time, ergsig2_hwrst); title('Photodiode2 Signal (ERG2):
HWR');
set(gca,'YLim',[0 1.5])
```





FIND ONSETS (DIFF(d_hwr==1) AND OFFSETS (diff(D_hwr==-1)).

```
    erg1_diff = diff(ergsig1_hwrst); % 10D of D_hwr
    erg1_diff = cat(2,erg1_diff,0);
    erg2_diff = diff(ergsig2_hwrst);
    erg2_diff = cat(2,erg2_diff,0);

    [pks_sig1,minima_sig1,locs_pks_sig1,locs_min_sig1]=
    CREx_peakfinder(erg1_diff); %Call of function to locate peaks.
    [pks_sig2,minima_sig2,locs_pks_sig2,locs_min_sig2]=
    CREx_peakfinder(erg2_diff);

    if locs_min_sig1(1)<locs_pks_sig1(1) % If offset precedes an
onset in time for photodiode 1
        onoffsets_sig1 =
[time(locs_pks_sig1);time(locs_min_sig1(2:end))];
        durs_sig1 = onoffsets_sig1(:,2) - onoffsets_sig1(:,1); %
Calculate the photodiode signal duration.
        onoffsets_ergol{1,ergcnt} = cat(2,onoffsets_sig1,durs_sig1); %
Concatenate onset/offset and signal duration information.
        onsets_all1 = nan(size(time));
        offsets_all1 = nan(size(time));
        onsets_all1(erg1_diff == pks_sig1(1)) = 0;
        offsets_all1(erg1_diff== minima_sig1(1)) = 0;
        onsets_ergol{1,ergcnt} = onsets_all1;
        offsets_ergol{1,ergcnt} = offsets_all1;

    else
        onoffsets_sig1 = [time(locs_pks_sig1);time(locs_min_sig1)];
        durs_sig1 = onoffsets_sig1(:,2) - onoffsets_sig1(:,1);
        onoffsets_ergol{1,ergcnt} = cat(2,onoffsets_sig1,durs_sig1);
        onsets_all1 = nan(size(time));
        offsets_all1 = nan(size(time));
        onsets_all1(erg1_diff == pks_sig1(1)) = 0;
        offsets_all1(erg1_diff== minima_sig1(1)) = 0;
        onsets_ergol{1,ergcnt} = onsets_all1;
        offsets_ergol{1,ergcnt} = offsets_all1;
    end

    if locs_min_sig2(1)<locs_pks_sig2(1) % If offset precedes an
onset in time for photodiode 2

        onoffsets_sig2 =
[time(locs_pks_sig2);time(locs_min_sig2(2:end))];
        durs_sig2 = onoffsets_sig2(:,2) - onoffsets_sig2(:,1); %
Calculate the photodiode signal duration.
        onoffsets_ergo2{1,ergcnt} = cat(2,onoffsets_sig2,durs_sig2);
        onsets_all2 = nan(size(time));
```

```

        offsets_all2 = nan(size(time));
        onsets_all2(erg2_diff == pks_sig2(1)) = 0;
        offsets_all2(erg2_diff== minima_sig2(1)) = 0;
        onsets_ergo2{1,ergcnt} = onsets_all2;
        offsets_ergo2{1,ergcnt} = offsets_all2;
    else

        onoffsets_sig2 = [time(locs_pks_sig2);time(locs_min_sig2)]';
        durs_sig2 = onoffsets_sig2(:,2) - onoffsets_sig2(:,1);
        onoffsets_ergo2{1,ergcnt} = cat(2,onoffsets_sig2,durs_sig2);
        onsets_all2 = nan(size(time));
        offsets_all2 = nan(size(time));
        onsets_all2(erg2_diff == pks_sig2(1)) = 0;
        offsets_all2(erg2_diff== minima_sig2(1)) = 0;
        onsets_ergo2{1,ergcnt} = onsets_all2;
        offsets_ergo2{1,ergcnt} = offsets_all2;
    end

```

PLOT THE DETRENDED AND FILTERED PHOTODIODE SIGNAL (RIGHT) AND SIGNAL CONVERTED TO STEP FUNCTION (LEFT).

```

figure('Name',fb_sessions{1,ergcnt},'NumberTitle','off');
subplot(2,2,1)
plot(time,ergsig1_hwrst);
hold on
plot(time,onsets_ergo1{1,ergcnt},'or','MarkerFaceColor','r')
hold on
plot(time,offsets_ergo1{1,ergcnt},'og','MarkerFaceColor','g');
set(gca,'YLim',[0 1.5])
title('Photodiode Signal as Step-function: onsets and offsets');
subplot(2,2,2)
plot(time,ergsig1Dfilt_corr)
hold on
plot(time,onsets_ergo1{1,ergcnt},'or','MarkerFaceColor','r');
hold on
plot(time,offsets_ergo1{1,ergcnt},'og','MarkerFaceColor','g');
title('Photodiode signal (detrend+filter): onsets (red), offsets (green)');

subplot(2,2,3)
plot(time,ergsig2_hwrst);
hold on
plot(time,onsets_ergo2{1,ergcnt},'or','MarkerFaceColor','r')
hold on
plot(time,offsets_ergo2{1,ergcnt},'og','MarkerFaceColor','g');
set(gca,'YLim',[0 1.5])
title('Photodiode Signal as Step-function: onsets and offsets');
subplot(2,2,4)
plot(time,ergsig2Dfilt_corr)
hold on


```

```


plot(time,onsets_ergo2{1,ergcnt},'or','MarkerFaceColor','r');
hold on
plot(time,offsets_ergo2{1,ergcnt},'og','MarkerFaceColor','g');
title('Photodiode signal (detrend+filter): onsets (red), offsets
(green)');

```


Photodiode Signal as Step-function: onsets (red), offsets (green)




Photodiode Signal as Step-function: onsets (red), offsets (green)




Photodiode Signal as Step-function: onsets (red), offsets (green) Photodiode signal (detrend+filter): onsets (red), offsets (green)




Photodiode Signal as Step-function: onsets (red), offsets (green) Photodiode signal (detrend+filter): onsets (red), offsets (green)




Photodiode Signal as Step-function: onsets (red), offsets (green) Photodiode signal (detrend+filter): onsets (red), offsets (green)



Photodiode Signal as Step-function: onsets (red), offsets (green) Photodiode signal (detrend+filter): onsets (red), offsets (green)





PLOT THE DETAIL OF THE PHOTODIODE SIGNALS TO SHOW DETECTED ONSETS AND OFFSETS.

```


indx = 100000:120000; % Indices to be included in plot.





figure('Name',fb_sessions{1,ergcnt},'NumberTitle','off');
subplot(1,2,1)
plot(time(indx),ergsig1_hwrst(indx));
hold on
plot(time(indx),onsets_ergo1{1,ergcnt}
(indx),'or','MarkerFaceColor','r')
hold on
plot(time(indx),offsets_ergo1{1,ergcnt}
(indx),'og','MarkerFaceColor','g')
set(gca,'YLim',[0 1.5])


subplot(1,2,2)
plot(time(indx),ergsig2_hwrst(indx));
hold on
plot(time(indx),onsets_ergo2{1,ergcnt}
(indx),'or','MarkerFaceColor','r')
hold on

```

```
plot(time(indxs),offsets_ergo2{1,ergcnt}  
(indxs),'og','MarkerFaceColor','g')  
set(gca,'YLim',[0 1.5])
```







end

CATEGORIZE THE FEEDBACK ONSET TIMES BASED ON DURATIONS

Based on the calculated durations, determine the feedback modality (verbal, gestural) and type (congruent, incongruent).

```
triginfo_file =
    fullfile(filesep, 'Users', 'bolger', 'Documents', 'work', 'Projects', 'Project-
    PEPs', 'Brain_IHM_Lists.xlsx'); %Path to excel file with trigger info.

% Initialize variables.
onsets_all = cell(1, size(ergsig1, 2));
fbtypes_all = cell(1, size(ergsig1, 2));
lats_all = cell(1, size(ergsig1, 2));
onset_data = cell(1, size(ergsig1, 2));
Allergo_onsets = cell(1, size(ergsig1, 2));

bloc_num = length(filenom2);
durs_erg1onset = cell(1, bloc_num);
durs_erg2onset = cell(1, bloc_num);
ideg1=cell(1, bloc_num);
ideg2=cell(1, bloc_num);
outlier_eg1=cell(1, bloc_num);
```

```

outlier_eg2=cell(1,bloc_num);
trialnumbers = cell(1,bloc_num);
fb_table = cell(1,bloc_num);

trignames = {'Gestual-cong' 'Verbal-cong' 'Gestual-incong' 'Verbal-
incong'}; % The 4 conditions.

for bcnt = 1:bloc_num

    % Read in information on Lists, corresponding videos and trigger-
types
    % from xlsx file.
    ix = strfind(filenum2{1,bcnt},'_');
    sheetnom = filenum2{1,bcnt}(ix(1)+1:ix(2)-1);
    XIn =
readtable(triginfo_file,'FileType','spreadsheet', 'sheet',sheetnom, 'ReadVariableNames',
'TreatAsEmpty','NA', 'Range','A1:F5');
    display(XIn); %print to screen a table showing the properties of
current list.
    fprintf('The current film is:\t%s\n', char(XIn{bcnt,1})); %print
to screen the current film

```

Warning: Variable names were modified to make them valid MATLAB identifiers. The original names are saved in the VariableDescriptions property.

XIn =

4x5 table

ERG2_Congruent	Film_Type		ERG1_Congruent	
	ERG1_Incongruent		ERG2_Incongruent	
Film1	'Curare_Human_Congruent'		200	160
	NaN	NaN		
Film2	'Perf_Agent_Congruent'		200	160
	NaN	NaN		
Film3	'Reg_Agent_Incongruent'		200	160
	400	320		
Film4	'Curare_Human_Incongruent'		200	160
	400	320		

The current film is: Curare_Human_Congruent

Warning: Variable names were modified to make them valid MATLAB identifiers. The original names are saved in the VariableDescriptions property.

XIn =

4x5 table

ERG2_Congruent	Film_Type	ERG1_Congruent	
	ERG1_Incongruent	ERG2_Incongruent	
Film1	'Curare_Human_Congruent'	200	160
	NaN	NaN	
Film2	'Perf_Agent_Congruent'	200	160
	NaN	NaN	
Film3	'Reg_Agent_Incongruent'	200	160
	400	320	
Film4	'Curare_Human_Incongruent'	200	160
	400	320	

The current film is: Perf_Agent_Congruent

Warning: Variable names were modified to make them valid MATLAB identifiers. The original names are saved in the VariableDescriptions property.

XIn =

4x5 table

ERG2_Congruent	Film_Type	ERG1_Congruent	
	ERG1_Incongruent	ERG2_Incongruent	
Film1	'Curare_Human_Congruent'	200	160
	NaN	NaN	
Film2	'Perf_Agent_Congruent'	200	160
	NaN	NaN	
Film3	'Reg_Agent_Incongruent'	200	160
	400	320	
Film4	'Curare_Human_Incongruent'	200	160
	400	320	

The current film is: Reg_Agent_Incongruent

Warning: Variable names were modified to make them valid MATLAB identifiers. The original names are saved in the VariableDescriptions property.

XIn =

4x5 table

ERG2_Congruent	Film_Type	ERG1_Congruent	
	ERG1_Incongruent	ERG2_Incongruent	

Film1	'Curare_Human_Congruent'	200	160
	NaN	NaN	
Film2	'Perf_Agent_Congruent'	200	160
	NaN	NaN	
Film3	'Reg_Agent_Incongruent'	200	160
	400	320	
Film4	'Curare_Human_Incongruent'	200	160
	400	320	

The current film is: Curare_Human_Incongruent

DETECT THE TRIGGERS BASED ON VIDEO-TYPE CORRESPONDING TO THE CURRENT BLOCK.

```

alltrigs = XIn{bcnt,2:5};
itrig = ~isnan(XIn{bcnt,2:5});
trigs_curr = alltrigs(itrig);
trignom_curr = trignames(itrig);

EG1durs = onoffsets_ergol{1,bcnt}(:,3); %durations of photodiode
signal for ergol
EG2durs = onoffsets_ergo2{1,bcnt}(:,3); %durations of photodiode
signal for ergo2
if length(trigs_curr)==2
    ideg1{1,bcnt} = EG1durs<=160 | EG1durs>=320;
    ideg2{1,bcnt} = EG2durs<=80 | EG2durs>=200;
elseif length(trigs_curr)==4
    ideg1{1,bcnt} = EG1durs<=160 | EG1durs>=460;
    ideg2{1,bcnt} = EG2durs<=80 | EG2durs>=400;
end
durs_erglonset{1,bcnt} = EG1durs;
durs_erg2onset{1,bcnt} = EG2durs;
outlier_eg1{1,bcnt} = EG1durs(ideg1{1,bcnt});
outlier_eg2{1,bcnt} = EG2durs(ideg2{1,bcnt});
% onoffsets_ergol{:,bcnt} = onoffsets_ergol{:,bcnt}
(~ideg1{1,bcnt});
% onoffsets_ergo2{:,bcnt} = onoffsets_ergo2{:,bcnt}
(~ideg2{1,bcnt});

%ERG01 (gestual) to begin with...
if length(trigs_curr)<=2
    trig_erglincong = 400;
    trig_erg2incong = 320;
else
    trig_erglincong = trigs_curr(3);
    trig_erg2incong = trigs_curr(4);
end

erglincong_diff = arrayfun(@(eg1) abs(eg1-trig_erglincong),
durs_erglonset{1,bcnt});

```

```

    erg1cong_diff = arrayfun(@(eg1) abs(eg1-trigs_curr(1)),
durs_erglonset{1,bcnt});
    [erg1min, ierg1] = min([erg1cong_diff, erg1incong_diff],[],2);
    % 1=congruent; 2=incongruent

    %ERGO2 (verbal)...
    erg2incong_diff = arrayfun(@(eg2) abs(eg2-trig_erg2incong),
durs_erg2onset{1,bcnt});
    erg2cong_diff = arrayfun(@(eg2) abs(eg2-trigs_curr(2)),
durs_erg2onset{1,bcnt});
    [erg2min, ierg2] = min([erg2cong_diff, erg2incong_diff], [],2);
    % 1=congruent; 2=incongruent

    % Determine the latencies for the ERGO1 and ERGO2 onsets
    time = Times_bloc{bcnt,1};
    erglonsets = onoffsets_ergo1{:,bcnt}(:,1);
    erglonsets_cong = erglonsets(ierg1==1); %ergo1 congruent onset
times
    erglonsets_incong = erglonsets(ierg1==2); %ergo1 incongruent onset
times

    erg2onsets = onoffsets_ergo2{:,bcnt}(:,1);
    erg2onsets_cong = erg2onsets(ierg2==1); %ergo2 congruent onset
times
    erg2onsets_incong = erg2onsets(ierg2==2); %ergo2 incongruent onset
times

    if isempty(erglonsets_incong) || isempty(erg2onsets_incong)

        dumerg1_cong = ones(length(erglonsets_cong),1);
        dumerg2_cong = ones(length(erg2onsets_cong),1).*2;
        Allergs = [erglonsets_cong; erg2onsets_cong];
        Alldums = [dumerg1_cong; dumerg2_cong];
        all_idegs = [ideg1{1,bcnt};ideg2{1,bcnt}];
        dum_erg = [Allergs, Alldums];
        [onsets_all{1,bcnt},idx] = sort(dum_erg(:,1),'ascend');
        idum_all = Alldums(idx,:); %1=vis_cong, 2=aud_cong
        idegs_all = all_idegs(idx,:);
    else

        dumerg1_cong = ones(length(erglonsets_cong),1);
        dumerg2_cong = ones(length(erg2onsets_cong),1).*2;
        dumerg1_incong = ones(length(erglonsets_incong),1).*3;
        dumerg2_incong = ones(length(erg2onsets_incong),1).*4;
        Allergs = [erglonsets_cong; erg2onsets_cong;
erglonsets_incong; erg2onsets_incong];
        Alldums = [dumerg1_cong; dumerg2_cong; dumerg1_incong;
dumerg2_incong];
        all_idegs = [ideg1{1,bcnt};ideg2{1,bcnt}];
        dum_erg = [Allergs, Alldums];
        [onsets_all{1,bcnt},idx] = sort(dum_erg(:,1),'ascend');
        idum_all = Alldums(idx,:); %1=vis_cong, 2=aud_cong,
3=vis_incong, 4=aud_incong

```

```

        idegs_all = all_idegs(idx,:);

end

    tdiff = arrayfun(@(tin) abs(tin-
onsets_all{1,bcnt})),time, 'UniformOutput', false);
    idx_t = cellfun(@(l) sum(l==0), tdiff);
    Allergo_onsets = time(idx_t==1);
    fbtype = cell(length(Allergo_onsets),1);
    lats = find(idx_t);

    %Assign trigger names to the onset times
    for i = 1:length(Allergo_onsets)

        if idegs_all(i)==1
            fbtype{i,1}=strcat(trignames{1,idum_all(i)},'-outlier');
        elseif ~ismember(trignames{1,idum_all(i)}, trignom_curr)
            fbtype{i,1} = strcat(trignames{1,idum_all(i)},'-outlier');
        else
            fbtype{i,1} = trignames{1,idum_all(i)};
        end
    end

end

fbtypes_all{1,bcnt} = fbtype;
lats_all{1,bcnt} = Allergo_onsets;

onset_data{1,bcnt}.video = string(fb_sessions{1,bcnt});
onset_data{1,bcnt}.onset_times = onsets_all{1,bcnt}';
onset_data{1,bcnt}.latencies = lats_all{1,bcnt}';
onset_data{1,bcnt}.types = fbtypes_all{1,bcnt}';

for cntr = 1:length(lats_all{1,bcnt})
    ALLEEG(bcnc).event(cntr).type = string(fbtypes_all{1,bcnt}
(cntr));
    ALLEEG(bcnc).event(cntr).latency = lats(cntr);
    ALLEEG(bcnc).event(cntr).urevent = cntr;
    ALLEEG(bcnc).urevent(cntr).type = string(fbtypes_all{1,bcnt}
(cntr));
    ALLEEG(bcnc).urevent(cntr).latency = lats(cntr);
end

[ALLEEG, ~, CURRENTSET] = pop_newset(ALLEEG, ALLEEG(bcnc),
CURRENTSET,'setname',char(fb_sessions{1,bcnt}),'gui','off');
[ALLEEG, EEG] = eeg_store(ALLEEG, ALLEEG(bcnc), CURRENTSET);
EEG = eeg_checkset( EEG );
EEG =
pop_saveset( EEG, 'filename',char(fb_sessions{1,bcnt}),'filepath',Savebase);
eeglab redraw

```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 5

Warning: converting all event types to strings
Saving dataset...

#5: S01_List1a_Film1

```
Filename: ...ts/s01/S01 List1a Film1.set
Channels per frame           74
Frames per epoch             563200
Epochs                       1
Events                       98
Sampling rate (Hz)           2048
Epoch start (sec)            0.000
Epoch end (sec)              275.000
Reference                     unknown
Channel locations             No (labels only)
ICA weights                   No
Dataset size (Mb)             171.4
```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 7
Warning: converting all event types to strings
Saving dataset...

#7: S01_List1a_Film2

```
Filename: ...ts/s01/S01 List1a Film2.set
Channels per frame           74
Frames per epoch             403456
Epochs                       1
Events                       86
Sampling rate (Hz)           2048
Epoch start (sec)            0.000
Epoch end (sec)              197.000
Reference                     unknown
Channel locations             No (labels only)
ICA weights                   No
Dataset size (Mb)             122.8
```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 9
Warning: converting all event types to strings
Saving dataset...

#9: S01_List1a_Film3

```
Filename: ...ts/s01/S01 List1a Film3.set
Channels per frame          74
Frames per epoch           399360
Epochs                     1
Events                      64
Sampling rate (Hz)         2048
Epoch start (sec)          0.000
Epoch end (sec)            195.000
Reference                   unknown
Channel locations           No (labels only)
ICA weights                 No
Dataset size (Mb)           121.5
```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 11
Warning: converting all event types to strings
Saving dataset...

#11: S01_List1a_Film4

```
Filename: ...ts/s01/S01 List1a Film4.set
Channels per frame          74
Frames per epoch           477184
Epochs                     1
Events                      85
Sampling rate (Hz)         2048
Epoch start (sec)          0.000
Epoch end (sec)            233.000
Reference                   unknown
Channel locations           No (labels only)
ICA weights                 No
Dataset size (Mb)           145.2
```

CALCULATE ITI AND DETERMINE FEEDBACKS FALLING INTO THE SAME TRIAL.

```
min_dist = 0.1;    %define the minimum distance between two
feedbacks. Below this value, the feedbacks are considered as being in
the same trial.
diff_lat = diff([ALLEEG(bcnt).event.latency]);
```

```

diffflat_sec = diff_lat./ALLEEG(bcnt).srate;
trl_dummy = zeros(length(diffflat_sec)+1,1);
trlcounter = 1;

for trlcnt = 1:length(ALLEEG(bcnt).event)

    if trlcnt==1
        t_end=trlcounter+1;
    end

    if t_end<=length(ALLEEG(bcnt).event)

        diff_curr = (ALLEEG(bcnt).event(t_end).latency-
ALLEEG(bcnt).event(trlcounter).latency)/ALLEEG(bcnt).srate;

        if diff_curr<=min_dist
            trl_dummy(trlcounter:t_end)=trlcnt;
            trlcounter = t_end+1;
            t_end=t_end+2;

        elseif diff_curr> min_dist
            trl_dummy(trlcounter)=trlcnt;
            trlcounter = trlcounter+1;
            t_end=t_end+1;
        end

        elseif t_end>length(ALLEEG(bcnt).event) && diff_curr<=min_dist

            trl_dummy(end) = trl_dummy(end-1);

        elseif t_end>length(ALLEEG(bcnt).event) && diff_curr> min_dist

            trl_dummy(end)=trl_dummy(end-1)+1;
        end

        trl_dummy(end)=trl_dummy(end-1)+1;

    end

    trialnumbers{1,bcnt}=trl_dummy;

```

ADD TRIAL NUMBER INFORMATION TO THE EVENT FIELD OF THE EEG STRUCTURE (EEGLAB based).

```

for icnt = 1:length([ALLEEG(bcnt).event])

    ALLEEG(bcnt).event(icnt).trialnum = trl_dummy(icnt);

end

```

```

[ALLEEG, ~, CURRENTSET] = pop_newset(ALLEEG, ALLEEG(bcnt),
CURRENTSET, 'setname', char(fb_sessions{1,bcnt}),'gui','off');
[ALLEEG, EEG] = eeg_store(ALLEEG, ALLEEG(bcnt), CURRENTSET);
EEG = eeg_checkset( EEG );
EEG =
pop_saveset( EEG, 'filename',char(fb_sessions{1,bcnt}),'filepath',Savebase);
eeglab redraw

```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 6
Warning: converting all event types to strings
Saving dataset...

#6: S01_List1a_Film1

```

Filename: ...ts/s01/S01 List1a Film1.set
Channels per frame           74
Frames per epoch             563200
Epochs                       1
Events                       98
Sampling rate (Hz)           2048
Epoch start (sec)            0.000
Epoch end (sec)              275.000
Reference                     unknown
Channel locations             No (labels only)
ICA weights                   No
Dataset size (Mb)             171.4

```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 8
Warning: converting all event types to strings
Saving dataset...

#8: S01_List1a_Film2

```
Filename: ...ts/s01/S01 List1a Film2.set
Channels per frame          74
Frames per epoch            403456
Epochs                      1
Events                       86
Sampling rate (Hz)          2048
Epoch start (sec)           0.000
Epoch end (sec)             197.000
Reference                    unknown
Channel locations            No (labels only)
ICA weights                  No
Dataset size (Mb)            122.8
```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 10
Warning: converting all event types to strings
Saving dataset...

#10: S01_List1a_Film3

```
Filename: ...ts/s01/S01 List1a Film3.set
Channels per frame          74
Frames per epoch            399360
Epochs                      1
Events                       64
Sampling rate (Hz)          2048
Epoch start (sec)           0.000
Epoch end (sec)             195.000
Reference                    unknown
Channel locations            No (labels only)
ICA weights                  No
Dataset size (Mb)            121.5
```

Warning: converting all event types to strings
Creating a new ALLEEG dataset 12
Warning: converting all event types to strings
Saving dataset...

#12: S01_List1a_Film4

```
Filename: ...ts/s01/S01 List1a Film4.set
Channels per frame          74
Frames per epoch            477184
Epochs                      1
Events                      85
Sampling rate (Hz)          2048
Epoch start (sec)           0.000
Epoch end (sec)             233.000
Reference                    unknown
Channel locations            No (labels only)
ICA weights                  No
Dataset size (Mb)            145.2
```

SAVE ALL THE INFORMATION INTO A MAT FILE AS A TABLE FOR EACH FILM

```
Trial_numbers = trl_dummy;
Feedbacks = {EEG.event.type}';
Times = [ALLEEG(bcnt).event.latency]'./ALLEEG(bcnt).srate;
ITI = cat(1,difflat_sec',nan);

fb_table{1,bcnt} = table(Times, Feedbacks, Trial_numbers, ITI);

save(fullfile(Savebase,'feedback_summary.mat'),'fb_table');

eeglab redraw
```

#6: S01_List1a_Film1

```
Filename: ...ts/s01/S01 List1a Film1.set
Channels per frame          74
Frames per epoch           563200
Epochs                     1
Events                     98
Sampling rate (Hz)         2048
Epoch start (sec)          0.000
Epoch end (sec)            275.000
Reference                   unknown
Channel locations           No (labels only)
ICA weights                 No
Dataset size (Mb)           171.4
```

#8: S01_List1a_Film2

```
Filename: ...ts/s01/S01 List1a Film2.set
Channels per frame          74
Frames per epoch           403456
Epochs                     1
Events                     86
Sampling rate (Hz)         2048
Epoch start (sec)          0.000
Epoch end (sec)            197.000
Reference                   unknown
Channel locations           No (labels only)
ICA weights                 No
Dataset size (Mb)           122.8
```

#10: S01_List1a_Film3

```
Filename: ...ts/s01/S01 List1a Film3.set
Channels per frame           74
Frames per epoch             399360
Epochs                       1
Events                       64
Sampling rate (Hz)           2048
Epoch start (sec)            0.000
Epoch end (sec)              195.000
Reference                     unknown
Channel locations             No (labels only)
ICA weights                   No
Dataset size (Mb)             121.5
```

#12: S01_List1a_Film4

```
Filename: ...ts/s01/S01 List1a Film4.set
Channels per frame           74
Frames per epoch             477184
Epochs                       1
Events                       85
Sampling rate (Hz)           2048
Epoch start (sec)            0.000
Epoch end (sec)              233.000
Reference                     unknown
Channel locations             No (labels only)
ICA weights                   No
Dataset size (Mb)             145.2
```

end

Published with MATLAB® R2018b